

Monivent® Neo Training

essential support for manikin based ventilation training

TECHNICAL SPECIFICATION*

Mask sizes S, M, L	40 mm, 52 mm, 66 mm	Size	Charging station: 267 x 80 x 65 mm Tablet: 169.5 x 240 x 7.5 mm Sensor module: 54 x 20 x 23 mm
Measured parameters	<ul style="list-style-type: none">• VTe• PIP• PEEP• Face mask leak• Ventilation rate	Wireless communication	Sensor module: Wireless data transfer to tablet Export data: WiFi
Sensor module operating time	1 hour		
Weight	Charging station: 475 g Tablet: 469 g Sensor module: 15 g		

* Technical specification is subject to change

For ordering or more information, please contact:

Manufacturer

Monivent AB
Regnbågsgatan 8B
417 55 Göteborg, Sweden

Tel: +46 31 395 48 40, E-post: info@monivent.se

monivent.se



Local distributor

Not yet available in all markets.
Please contact Monivent or local representative
for more information.



Monivent **Neo Training** offers a unique training solution for maintaining and improving manual ventilation skills.




Monivent® Neo Training

The ability to deliver effective ventilations represents a critical factor when managing newborns requiring ventilatory support.

As the need to apply this skill can be relatively infrequent, there is a clear need for regular and high quality ventilation training.

Monivent Neo Training is easy to set up, portable and provides a unique simulation setting. Practice at your own convenience, when and where you have time/ want. Every training session is logged and stored in the system to allow tracking of performance over time and compare results to previous sessions.



Monivent Neo training includes:

- 6 face masks (2 Small, 2 Medium, 2 Large)
- 3 sensor modules
- Charging station
- Application on tablet for feedback display
- Carrying case



Benefits

- + Simple set up by connecting the Monivent Neo Training face mask to conventional bag valve mask or T-piece ventilation device.
- + Wireless data transfer provides greater flexibility and freedom to operate by eliminating connecting cables or tubing.
- + Continuous real time feedback enhances and speeds up learning.
- + Feedback on face mask leak helps identify need to adjust face mask position.
- + PIP in combination with VT_e helps identify potential airway obstruction and provides a better understanding of lung compliance.
- + Personal training log showing performance over time, increases motivation to practice and improves ventilation skills.
- + Monivent Neo Training is compatible with all bag or T-piece resuscitators utilizing standard ISO 5356-1 15 mm female connection.
- + For accurate measurements and high-quality training a leak-free manikin with realistic lung compliance is required.



Application on tablet for feedback display

- Color-coded visual display of accurate ventilation with numerical indication or values of PIP, PEEP, leak and ventilation rate.



Lightweight wireless sensor module

- Wireless transfer of data allows greater freedom and flexibility to practice (no need to connect cables for data transfer) (placing monitor wherever wanted)
- Integrated color LED gives guidance on VT_e on target (green light) or outside pre-set target interval (orange/ red light)



Face masks

- Face masks in three different sizes: S, M, L
- Compatible with all bag or T-piece resuscitators utilizing standard connection



Charging station

- LED light indicates battery status
- Sensor modules and masks can be stored in charging station between use